

Vioron®

The versatile trypan blue dye
for the anterior segment



Composition and properties

1 ml of Vioron® contains:

- 0.6 mg Trypan Blue
- Disodium phosphate ($\text{Na}_2\text{HPO}_4 \times 2 \text{ H}_2\text{O}$)
- Monosodium phosphate ($\text{NaH}_2\text{PO}_4 \times 2 \text{ H}_2\text{O}$)
- Sodium chloride (NaCl)
- Water for injection

Density: 1.00 g/cm³

Brilliant visualization
of the anterior lens capsule

Excellent distinction of the
capsulorhexis margin

Quick and easy application
(Ready-to-use)

Approved for DMEK

Packaging units

G-81002 Vioron® Syringe
0.5 ml syringe, 5 pcs. per box, sterile

G-81001 Vioron® Vial
0.5 ml vial, 5 pcs. per box, sterile

Fields of application

Voron® was developed for ophthalmic surgical procedures in the anterior segment of the eye such as cataract operations or keratoplasties. Staining the anterior lens capsule makes it more visible, thus facilitating capsulorhexis and minimizing the risk of tearing. Furthermore, Voron® facilitates the preparation and transfer of the donor cornea in the case of lamellar corneal transplantsations and the removal of the diseased Descemet's membrane in case of DMEK and DS(A)EK.



References

Keratoplasty

49. Melles, G.R.J., T. San Ong, B. Ververs and J. van der Wees. 2006. "Descemet Membrane Endothelial Keratoplasty (DMEK)." *Cornea* 25(8): 987-990. 50. Ham L., I. Dapena, C. van Luijk, J. van der Wees and G. R. J. Melles. 2009. "Descemet Membrane Endothelial Keratoplasty (DMEK) for Fuchs endothelial dystrophy: review of the first 50 consecutive cases." *Eye* 23(10): 1990-1998. 51. Cursiefen, C. and F. E. Kruse. 2010. "DMEK: Posteriore lamellare Keratoplastik-technik." *Der Ophthalmologe* 107(4): 370-376. 52. Droutsas, K., L. Ham, I. Dapena, G. Geerling, S. Oellerich and G.R.J. Melles. 2010. "Visus nach Descemet-Membran Endothelkeratoplastik (DMEK): Ergebnisse der ersten 100 Eingriffe bei Fuchs'scher Endotheldystrophie." *Klinische Monatsblätter für Augenheilkunde* 227(06): 467-477. 53. Dapena, I., K. Moutsouris, K. Droutsas, L. Ham, K. van Dijk and G.R.J. Melles. 2011. "Standardized "To-Touch" Technique for Descemet Membrane Endothelial Keratoplasty." *Archives of Ophthalmology* 129(1): 88-94. 54. Dirisamer, M., L. Ham, I. Dapena, K. Moutsouris, K. Droutsas, K. van Dijk, L. E. Frank, S. Oellerich and G.R.J. Melles. 2011. "Efficiency of Descemet Membrane Endothelial Keratoplasty: Clinical Outcome of 200 Consecutive Cases after a Learning Curve of 25 Cases." *Archives of Ophthalmology* 129(11): 1435-1443. 55. Heindl, L. M., S. Riss, B. O. Bachmann, K. Laaser, F. E. Kruse and C. Cursiefen. 2011. "Split cornea transplantation for 2 recipients. A new strategy to reduce corneal tissue cost and shortage." *Ophthalmology* 118(2): 294-301. 56. Güell, J. L., M. Morral, O. Gris, D. Elies and F. Manero. 2013. "Bimanual technique for insertion and positioning of Endothelial-Descemet Membrane graft in Descemet Membrane Endothelial Keratoplasty." *Cornea* 32(12): 1522-1526. 57. Ciechanowski, P. P., K. Droutsas, L. Baydoun, M. Dirisamer, S. Oellerich, G.R.J. Melles. 2014. "Standardisierte Descemet-Membran-Endothelkeratoplastik (DMEK). Technik und aktuelle Ergebnisse." *Der Ophthalmologe* 111(11): 1041-1049. 58. Baydoun, L., K. van Dijk, I. Dapena, F. U. Musa, V. Liarakos, L. Ham and G.R.J. Melles. 2015. "Repeat Descemet Membrane Endothelial Keratoplasty after complicated primary Descemet Membrane Endothelial Keratoplasty." *Ophthalmology* 122(1): 8-16.

Capsulorhexis

59. Melles, G.R.J., P. W. T. de Waard, J. H. Pameijer and W. H. Beekhuis. 1999. "Trypan blue capsule staining to visualize the capsulorhexis in cataract surgery." *Journal of Cataract and Refractive Surgery* 25(1):7-9. 60. Melles, G.R.J., P. W. T. de Waard, J. H. Pameijer and W. H. Beekhuis. 1999. "Farbung der Linsenkapsel mit Trypanblau zur Visualisierung der Kapsulorhexis bei Maturkataraktchirurgie." *Klinische Monatsblätter für Augenheilkunde* 215(12): 342-344. 61. Werner, L., S. K. Pandey, M. Escobar-Gomez, D. S. Hoddinott and D. J. Apple. 2000. "Dye-enhanced cataract surgery. Part 2: Learning critical steps of phacoemulsification1." *Journal of Cataract and Refractive Surgery* 26(7): 1060-1065. 62. De Waard, P.W.T., C. J. Budo and G.R.J. Melles. 2002. "Trypan blue capsular staining to "find" the leading edge of a "lost" capsulorhexis." *American Journal of Ophthalmology* 134(2): 271-272. 63. Jacob, S., A. Agarwal, A. Agarwal, S. Agarwal, S. Chowdhary, R. Chowdhary and A. A. Bagmar. 2002. "Trypan blue as an adjunct for safe phacoemulsification in eyes with white cataract." *Journal of Cataract and Refractive Surgery* 28(10): 1819-1825. 64. Ozkiris, A., O. Arslan, E. Cicik, N. Köylüoglu and C. Evereklioglu. 2003. "Open-sky capsulorhexis in triple procedure: with or without trypan blue?" *European Journal of Ophthalmology* 13(9-10): 764-769. 65. Nardi, M., U. Benelli, M. Figus and M. P. Bartolomet. 2004. "Open-sky capsulorhexis in triple procedure: with or without trypan blue?" *European Journal of Ophthalmology* 14(4): 358-358. 66. Singh, A. J., U. A. Sarodia, L. Brown, R. Jagjivan and R. Sampath. 2003. "A histological analysis of lens capsules stained with trypan blue for capsulorhexis in phacoemulsification cataract surgery." *Eye* 17(5): 567-570. 67. Xiao, Y., Y. H. Wang, Z. Y. Fu and H. Hong. 2004. "Staining the anterior capsule with indocyanine green or trypan blue for capsulorhexis in eyes with white cataract." *International Ophthalmology* 25(5-6): 273-276. 68. Chung, C. F., C. C. Liang, J. S. Lai, E. S. Lo and D. S. Lam. 2005. "Safety of trypan blue 1 % and indocyanine green 0.5 % in assisting visualization of anterior capsule during phacoemulsification in mature cataract." *Journal of Cataract and Refractive Surgery* 31(5): 938-942. 69. Wong, V. W. Y., T. Y. Y. Lai, G. K. Y. Lee, P. T. H. Lam and D. S. C. Lam. 2006. "A prospective study on trypan blue capsule staining under air vs under viscoelastic." *Eye* 20(7): 820-825. 70. Chéour, M., F. B. Brahim, A. Zarrad, N. Khémiri, K. Mgħaieħ and A. Kraiem. 2007. "Phacoemulsification des cataractes blanches en utilisant le bleu trypan." *Journal Français d' Ophthalmologie* 30(9): 914-917. 71. Ziakas, N. G., K. Boboridis, E. Nakos, D. Mikropoulos, V. Margaritis and A.G.P. Konstas. 2009. "Does the use of trypan blue during phacoemulsification affect the intraocular pressure?" *Canadian Journal of Ophthalmology* 44(3): 293-296. 72. Smith, E. F., R. U. Desai, A. Schrier, B. Enriquez and B. K. Purewal. 2010. "Trypan blue capsulorhexis." *Ophthalmology* 117(7): 1462-1462.e1. 73. Jaber, R., L. Werner, S. Fuller, S. C. Kavoussi, S. McIntyre, M. Burrow and N. Marmalis. 2012. "Comparison of capsulorhexis resistance to tearing with and without trypan blue dye using a mechanized tensile strength model." *Journal of Cataract and Refractive Surgery* 38(3): 507-512. 74. Giannaria, D., M. Giannotti, A. Scopelliti, G. Pellegrini and B. Giannotti. 2013. "Under-air staining of the anterior capsule using trypan blue with a 30 g needle." *Clinical Ophthalmology* (Auckland, NZ), 7: 233-235.



Distributed by:

Medix21 Surgical

29 Jarden Mile Ngauranga 6035

p 0800 633 4921 / 021 734 762

e camille@medix21.co.nz

w medix21surgical.co.nz