



MADE IN JAPAN


FEATHER BIOPSY PUNCH

Feather Biopsy Punch offers excellent sharpness by ultra-precision sharpening technology which originates in surgical blade manufacturing.




Packaged with safety cover

Anti-roll hexagon design









Completely-sealed package by special film

To open, peel the package off by pulling apart the two flaps as indicated.



- Sterilized by Gamma radiation
- Single use only

Product No.	Inner Diameter	Package
FBP-15	1.5mm 	20pcs per Box
FBP-30	3mm 	
FBP-40	4mm 	
FBP-50	5mm 	
FBP-60	6mm 	
Actual size 		



FEATHER SAFETY RAZOR CO., LTD.
OVERSEAS TRADE DIVISION

3-70, OHYODO MINAMI 3-CHOME, KITA-KU, OSAKA 531-0075, JAPAN
PHONE: +81-6-6458-1638 FAX: +81-6-6458-1611
URL <http://www.feather.co.jp/> E-mail overseas@feather.co.jp



Distributed by:
Medix21 Surgical
29 Jarden Mile Ngauranga 6035
p 0800 633 4921 / 021 734 762
e camille@medix21.co.nz
w medix21surgical.co.nz