



Titanium Ophthalmic Instruments



[www.surgitrac.com](http://www.surgitrac.com)



## About Surgitrac Instruments

We are Surgitrac Instruments we manufacture specialist surgical instruments. Our instruments, sets and procedure packs for use in Ophthalmology are CE certified and comply with international standards.


We also offer single-use and reusable, stainless steel instruments.

## Titanium reusable microsurgical instruments

Surgitrac Instruments premium range of reusable precision quality Titanium microsurgical instruments feature:

- Anti-glare blue finish.
- Reusable – can be sterilised by all known techniques.
- Lightweight, Titanium is a third of the weight of stainless steel.
- Non-toxic.
- Non-corrosive.
- Anti-magnetic.
- Tungsten carbide coated tips for extra grip.

This brochure features frequently used Titanium microsurgical instruments for use in Ophthalmology. If you require an ophthalmic instrument that does not feature in this brochure, please contact us.



Titanium



Tying Forceps

**Castroviejo Tying Forceps**




**Code**

1401031  
1401032

**Shaft**

Straight  
Angled

5.0 mm long smooth jaws for 8/0 to 11/0 suture. Shafts available straight or angled 45°, Castroviejo handle. Length 108 mm.



**Curved Tying Forceps**




**Code**

1401018  
1401026

**Length**

85 mm  
115 mm

Curved shafts with 6 mm tying platforms. Ideal for 10/0 or 11/0 suture.



**Girard Tying Forceps**




**Code**

1401041  
1401042

**Shafts**

Straight  
Curved

4.5 mm tying platforms. Round knurled handle with guide pin. Length 105 mm.



**McPherson Tying Forceps**



**Code**

1401013  
1401014  
1401015  
1401016


**Shaft**

Standard  
Ultrafine  
Ultrafine  
Ultrafine

**Tips**

12 mm  
7 mm  
10 mm  
12 mm

Angled shafts with 4.5 mm tying platforms. Length 85 mm.



**McPherson Tying Forceps**




**Code**

1401023  
1401024

**Length**

85 mm  
115 mm

Angled 45° shafts with 4.5 mm tying platforms. Length 115 mm.




### Straight Tying Forceps



**Code**  
1401011  
1401012

**Length**  
85 mm  
115 mm




Straight shafts with 4.5 mm tying platforms.

### Straight Tying Forceps



**Code**  
1401021  
1401022

**Length**  
85 mm  
115 mm




Straight shafts with 4.5 mm tying platforms.

### Toothed Forceps

### Bishop-Harmon Tissue Forceps



**Code**  
1401156




1 x 2 teeth, delicate. Length 85 mm.

### Bonn Toothed Forceps



**Code**  
1401160




1 x 2 teeth, 0.12 mm. Length 85 mm.

### Colibri Conjunctiva Toothed Forceps



**Code**  
1401117



1 x 2 teeth, 0.12 mm, with 5.0 mm tying platforms. Length 85 mm.

Titanium



### Colibri Toothed Forceps




#### Code

1401113  
1401123

#### Length

85 mm  
115 mm




1 x 2 teeth, 0.12 mm, with tying platforms.

### Colibri Toothed Forceps



#### Code

1401153



1 x 2 teeth, 0.12 mm, with tying platforms. Length 85 mm.

### Harms-Colibri Toothed Forceps




#### Code

1401151  
1401161

#### Teeth

0.12 mm  
0.30 mm




1 x 2 teeth, with tying platforms. Length 75 mm.

### Kelman-McPherson Toothed Forceps



#### Code

1401114




1 x 2 teeth, 0.12 mm, with 5.0 mm tying platforms and 10 mm angled shafts. Length 85 mm.

### Polack Corneal Suturing Forceps



#### Code

1401372




1 mm spread, very delicate 1 x 2 teeth, 0.12 mm. Length 115 mm.



### Polack Corneal Suturing Forceps

**Code**  
1401381





Colibri style. Two 2.75 mm long tips, separated 1 mm, with 1 x 2 teeth. Length 75 mm.

### Straight Toothed Forceps

**Code**  
1401111  
1401112

**Teeth**  
0.12 mm  
0.10 mm




1 x 2 teeth, with 5.0 mm tying platforms. Length 85 mm.

### Straight Toothed Forceps

**Code**  
1401121  
1401122


**Teeth**  
0.12 mm  
0.10 mm



1 x 2 teeth, with 5.0 mm tying platforms. Length 115 mm.

### Straight Toothed Forceps

**Code**  
1401141




1 x 2 teeth, 0.12 mm with tying platforms.

### Curved Toothed Forceps

**Code**  
1401152  
1401154

**Length**  
85 mm  
115 mm



1 x 2 teeth, 0.12 mm, with 5.0 mm tying platforms. Length 105 mm.

Titanium



### Tissue Forceps




1 x 2 teeth, delicate. Length 100 mm.

#### Code

1401159  
1401162

#### Shafts

Straight  
Curved



### Toothed Forceps




1 x 2 teeth, with 5.0 mm tying platforms. Length 108 mm.

#### Code

1401131  
1401132

#### Teeth

0.12 mm  
0.40 mm



### Notched Forceps


### Angled Notched Forceps



0.25 mm notched, with tying platforms. Length 85 mm.

#### Code

1401063




### Curved Notched Forceps



0.25 mm notched, with tying platforms. Length 85 mm.

#### Code

1401065




### Hofmann-Polack Corneal Suturing Forceps



1 mm spread, 0.25 mm notched. Length 115 mm.

#### Code

1401371





**Pierse Corneal Forceps**




0.25 mm notched, very delicate. Length 115 mm.

**Code**

1401081  
1401082

**Shafts**

Straight  
Curved




**Pierse-Colibri Corneal Forceps**



0.25 mm notched, with tying platforms. Length 75 mm.

**Code**

1401101



**Straight Notched Forceps**




0.25 mm notched, with tying platforms. Length 85 mm.

**Code**

1401061  
1401062

**Tips**

Standard  
Ultrafine



**Straight Notched Forceps**




0.25 mm notched. Length 115 mm.

**Code**

1401071  
1401072

**Shafts**

Straight  
Curved



Titanium




Capsulorhexis Forceps

**Bloomberg Pierser Capsulorhexis Forceps**

**Code**  
1401222



12 mm 45° angled curved shafts with an iris stop. With guide pin. Length 105 mm.




**Capsulorhexis Forceps**

**Code**  
1401225



45° angled 12 mm straight shafts and Iris stop. Length 120 mm.




**Capsulorhexis Forceps**

**Code**  
1401226



45° angled 12 mm curved shafts and Iris stop. Length 120 mm.




**Corydon Capsulorhexis Forceps**

**Code**  
1401223



Round handle. Fine 45° angled tips with curved 12 mm shafts and iris stop. Length 105 mm.




**Corydon Style Capsulorhexis Forceps**

<b>Code</b>	<b>Length</b>
1401176	85 mm
1401186	115 mm



Fine 45° angled sharp tips with curved 12 mm shafts.




**Lehner-Utrata Style Capsulorhexis Forceps**

**Code**  
1401221



Fine angled sharp tips with curved 12 mm shafts. With guide pin. Length 105 mm.




**O'Gawa-Utrata Capsulorhexis Forceps**

**Code**  
1401224



Fine 45° angled tips with shafts 12 mm with an iris stop Round handle. With guide pin. Length 105 mm.





**Utrata Style Capsulorhexis Forceps**

**Code**  
1401173  
1401174

**Shafts**  
Straight  
Curved



Fine angled sharp tips, 12 mm 45° angled shafts. Length 85 mm.




**Utrata Style Capsulorhexis Forceps**

**Code**  
1401227



Round handled. 45° angled tips with slim 10 mm shafts. Length 105 mm.



Titanium




IOL Forceps

**Lens Forceps**

Code	Length
1401261	85 mm
1401271	115 mm



Angled 45° shafts. With 7.5 mm long smooth jaws.




**Kelman-McPherson Lens Forceps**

Code
1401301



7.5 mm long smooth delicate jaws. Length 108 mm.




**Kratz Style Lens Forceps**

Code
1401283



Curved narrow smooth delicate jaws. Length 120 mm.




Chalazion Forceps

**Chalazion Forceps**

Code	Plate
1401405	Oval 12 x 14 mm
1401406	Round 8 mm
1401407	Round 15 mm
1401408	Oval 17 x 9 mm



Length 98 mm.




**Baby Jones Towel Clamps**

**Code**  
1401417



Length 65 mm.




**Cilia Forceps**

**Code**  
1401412



Length 72 mm.




**Hartmann Curved Hemostatic Mosquito Forceps**

**Code**  
1401415



Curved 20 mm serrated jaws. Length 95 mm.




**Hartmann Straight Hemostatic Mosquito Forceps**

**Code**  
1401414



Straight 20 mm serrated jaws, ring handle with ratchet lock. Length 95 mm.




**Scleral Plug Forceps**

**Code**  
1401416



Length 118 mm.




Titanium




Jewellers Forceps

**Jeweller Style Forceps**

**Code**  
1401421



0.1 x 0.06 mm tips. Length 110 mm.




**Angled Jeweller Style Forceps**

**Code**  
1401422




0.05 x 0.01 mm tips, 90° angled tips. Length 110 mm.




**Angled Jeweller Style Forceps (45 degree)**

**Code**  
1401423




0.05 x 0.01 mm tips bent at 45° angled tip. Length 110 mm.




**Curved Jeweller Style Forceps**

**Code**  
1401424



0.07 x 0.03 mm tips. Length 110 mm.






**Barraquer Solid Speculum**




**Code**  
1405008

With 14.5 mm long blades. Length 52 mm.




**Eye Speculum**



**Code**  
1405014

Swivel open blades to accommodate the variability of patient's anatomical features. Length 65 mm.




**Kratz Barraquer Titanium Wire Speculum**



**Code**  
1405007

With 13.5 mm open blades. Length 52 mm.




**Lieberman Speculum, K-wire**



**Code**  
1405005

Adjustable wire speculum with 14.5 mm open blades. With shaft angled 45°. Length 65 mm.




**Lieberman Speculum, V-wire**



**Code**  
1405002


The V-shaped wire blades are designed to help the surgeon retain the drape under the eye lids. Length 65 mm.



Titanium




### Naqahara Titanium Wire Speculum



With 13.5 mm blades. Length 52 mm.

**Code**  
1405009




### Reversible Lid Speculum



Length 65 mm.

**Code**  
1405006




### Shapira Speculum



Length 65 mm.

**Code**  
1405012




### Slade LASIK Titanium Wire Speculum




Length 65 mm.

**Code**  
1405003




### Titanium Wire Speculum



With 13.5 mm long blades. Length 52 mm.

**Code**  
1405001




**Titanium Wire Speculum**



Length 65 mm.

**Code**  
1405015




Titanium

**T** Nucleus Rotators & Choppers

**Surgitrac**  
Instruments


### Akahoshi Nucleus Stainer/manipulator

**Code**  
1404140




### Claw Phaco Chopper

**Code**  
1404130



### Drysdale Nucleus Manipulator


**Code**  
1404113



Paddle shaped tip.

### Haefliger Phaco Cleaver


**Code**  
1404139



Chopper and manipulator in one.


### Helveston Retractors

<b>Code</b>	<b>Tips</b>
1404158	7 mm
1404159	9 mm
1404160	11 mm




**Helveston Scleral Marking Ruler**

**Code**  
1404161




**Intracapsular Manipulator**

**Code**  
1404120




**Intracapsular Manipulator**

**Code**  
1404151



**Jaffe-Bechert Lens Nucleus Rotator**


**Code**  
1404105



Blunt forked tip.

**Lester IOL Manipulator**

<b>Code</b>	<b>Shaft</b>
1404106	Angled and curved
1404107	Straight and curved



IOL manipulator for lens with or without holes.

Titanium


**T** Nucleus Rotators & Choppers



### Manipulator



**Code**  
1404135




Angled 45° shaft, 10 mm from tip to bend.

### Nichamin Phaco Chopper



**Code**  
1404123  
1404124

**Blade**  
Straight  
Curved




Ambidextrous nucleus chopper, tip 1.1 mm long. Rounded and highly polished distal tip.

### Phaco Chopper



**Code**  
1404121  
1404122

**Shaft**  
Straight  
Angled




### Phaco Chopper



**Code**  
1404125  
1404125L

**Hand**  
Right  
Left




Phaco side cutter chopper, rounded and highly polished distal tip.

### Phaco Chopper



**Code**  
1404126  
1404126L

**Hand**  
Right  
Left



Phaco side cutter chopper, rounded and highly polished distal tip.




**Phaco Chopper**



**Code**  
 1404136  
 1404137

**Blade**  
 Straight  
 Curved




Ambidextrous nucleus chopper, tip 1.8 mm long Rounded and highly polished distal tip.

**Phaco Chopper**



**Code**  
 1404138  
 1404138L

**Hand**  
 Right  
 Left




Phaco side cutter chopper, rounded and highly polished distal tip. Tip length 1.8 mm.

**Phaco Chopper**



**Code**  
 1404141




Special inner blade is offset 45, rounded distal tip is highly polished, unique elongated distal segment.

**Phaco Chopper**



**Code**  
 1404142  
 1404142L

**Hand**  
 Right  
 Left




Rounded and highly polished distal tip. Tip length 1.8 mm.

**Phaco Chopper**



**Code**  
 1404302



Paddle manipulator / 60 Phaco chopper.


Titanium

**T** Nucleus  
Rotators &  
Choppers

**Surgitrac**  
Instruments

### Phaco Chopper

**Code**  
1404303




Double ended for vertical and horizontal phaco chopping.

### Repositor

**Code**  
1404114  
1404115


**Tips**  
0.7 mm  
0.5 mm



Synechia spatula. Semi-sharp tips.

### Schocket Double Ended Scleral Depressor

**Code**  
1404304



Double-end depressor and marker, 6 mm wide bar and 2.5 mm tear drop shape tips.

**Helveston Finder Hook**




Lightly angled cone shaped tip.

**Code**

1404155  
1404156

**Tips**

Small  
Large



**Helveston Muscle Hook**




Delicate muscle hook with blunt spatulated tip.

**Code**

1404152  
1404153  
1404154

**Tips**

8 mm  
10 mm  
12 mm




**Helveston Teaser Hook**



6 mm hook.

**Code**  
1404157




**Kuglen Iris Hook**



Iris hook and lens manipulator.

**Code**  
1404101  
1404102

**Shaft**  
Angled  
Straight




**Retinal Detachment Hook**



With 1.5 mm oval hole. Flattened 13 mm.

**Code**  
1404119



Titanium




### Sinsky Lens Manipulating Hook

**Code**  
1404110



0.15 mm dia. blunt tip, angled to right.



### Barraquer Needle Holder




**Code**

1402238  
1402245

**Lock**

Without  
With



Very delicate curved jaws, sharp tips. Length 115 mm.

### Barraquer Needle Holder-Curved




**Code**

1402101  
1402106

**Lock**

Without  
With



Delicate, 9 mm gently curved jaws, sharp tips. Length 115 mm.

### Bechart-Sinsky Needle Holder




**Code**

1402235  
1402243

**Lock**

Without  
With



Delicate, 9 mm gently curved jaws, sharp tips. Length 115 mm.

### Castroviejo Needle Holder




**Code**

1402191  
1402192

**Lock**

Without  
With



Gently curved delicate jaws, blunt tips. Length 140 mm.

### Castroviejo Needle Holder




**Code**

1402193  
1402195

**Lock**

With  
Without



Standard, 11 mm curved jaws, blunt tips. Length 140 mm.

Titanium

**Ti** Needle Holders



### Castroviejo Needle Holder




Delicate jaws. Straight blunt tips. Length 140 mm.

#### Code

1402204  
1402205

#### Lock

Without  
With



### Castroviejo Needle Holder




Standard, 11 mm jaws, blunt tips. Length 140 mm.

#### Code

1402203  
1402206

#### Lock

Without  
With



### Castroviejo Needle Holder




Fine 9 mm curved jaws, wide serrated handle, sharp tips. Length 125 mm.

#### Code

1402211  
1402212

#### Lock

Without  
With



### Castroviejo Needle Holder




Fine 9 mm straight jaws, wide serrated handle, sharp tips. Length 125 mm.

#### Code

1402213  
1402214

#### Lock

Without  
With



### Curved Needle Holder




Curved, 11 mm medium jaws, sharp tips. Length 115 mm.

#### Code

1402011  
1402016

#### Lock

Without  
With





### Curved Needle Holder




Curved, 6 mm delicate jaws. Length 110 mm.

#### Code

1402233  
1402241

#### Lock

Without  
With



### Fine Curved Needle Holder




Curved, 9 mm fine jaws, sharp tips. Length 115 mm.

#### Code

1402237  
1402244

#### Lock

Without  
With




### Side Curved Needle Holder



Blunt tips. Length 115 mm.

#### Code

1402236



### Straight Needle Holder




Straight, 6 mm delicate jaws. Length 110 mm.

#### Code

1402234  
1402242

#### Lock

Without  
With



### Troutman Needle Holder-Curved




10 mm gently curved jaws, blunt tips. Length 115 mm.

#### Code

1402231  
1402232

#### Lock

Without  
With




Titanium



### Fine Fixation Ring


**Code**  
1404149



Swivel head, 13 mm/9 mm cut out. With blunt teeth on both sides.

### Fine/Thornton Fixation Ring



<b>Code</b>	<b>Cut out</b>
1404147	13 mm/11 mm
1404148	16 mm/14 mm



Swivel head.

### Fixation Forceps


<b>Code</b>	<b>Spread</b>
1401373	13.0 mm
1401375	3.0 mm



1 into 2, 0.12 mm teeth at each tip. Length 115 mm.

### Helveston Fixation Forceps


**Code**  
1401382



1 x 2 0.5 mm teeth. Length 118 mm.

### Hofman-Thornton Fixation Ring

**Code**  
1404146




Swivel head, 16 mm diameter.

**Buratto Flap Protector**

**Code**  
1405031




Ideal for retracting the flap for maximum exposure of the bed and for protecting nasally hinged corneal flap during the laser ablation.




**Corneal Epithelium Spatula**

**Code**  
1405027



**Flap Elevator**

**Code**  
1405022




**Genisi Speculum LASIK**

**Code**  
1405034



Adjustable, strong with ideal curvature for accommodating microkeratome suction ring.  
Length 65 mm.




**Katzen LASIK Flap Unzipper**

**Code**  
1405046



Combines a modified Sinsky hook and a bi-convex spatula with bevelled notches for LASIK enhancements.




Titanium



### LASIK Flap Manipulator



**Code**  
1405045




Features a delicate Olive-shaped tip.

### Liberman Solid Speculum



**Code**  
1405029




Length 65 mm.

### Machat Adjustable Speculum



**Code**  
1405030




Length 65 mm.

### Steinert LASIK Speculum LASIK



**Code**  
1405033



Length 65 mm.

**Blade Breaker**

**Code**  
1410001



**Blade Handle**

**Code**  
1410002




**Braunstein Fixed Caliper**

<b>Code</b>	<b>Tips</b>
1410003	3.5/4.0 mm
1410004	3.0/3.5 mm



Fine tips mark sclerotomy sites for vitreous surgery




**Castroviejo Caliper**

**Code**  
1410005




Caliper measures from 0.0 mm to 20 mm in 0.5 mm increments with scale reading to one side.



**Christen Punctum Punch**

**Code**  
1410014




Titanium



### Descemet Membrane Punch


**Code**  
1410012



1.0 mm diameter head, 0.75 mm deep bite.

### Holth Corneo scleral Punch

**Code**  
1410013



1.0 mm wide by 7 mm long punch.

### Irrigating Vectis (3x8)


**Code**  
1410302



Cannula irrigating vectis, dual fluid flow through to front opening ports.

### Luntz-Dodick Trabulectomy Punch


<b>Code</b>	<b>Head Dia.</b>
1410011	1.0 mm
1410015	1.5 mm
1410016	1.2 mm



0.5 mm deep bite.

### Trabeculotome set

**Code**  
1404134




9 mm long pointed with 3 mm spread pair consists of left and right.

**Trabeculotome set**

**Code**  
1404134C



8 mm long pointed with 3 mm spread pair consists of left and right.




**Urrets-Zavalía sclera marker**

**Code**  
1404162



Pressure and mark the surface of sclera, blunt tips.





Reusable



Scissors



Please note: reusable scissors are supplied in stainless steel


Stainless Steel Scissors

### Westcott Scissors

Code  
UA11



Curved scissors.




### Ong Type Scissors

Code  
RUC80



Sharp pointed tips, angled blades.




### Vannas Scissors

Code  
RUC81



Straight scissors.




### Castroviejo Corneal Section Scissors (Right)

Code  
RUC82



13 mm from pivot to tip, right, 105 mm long.




### Castroviejo Corneal Section Scissors (Left)

Code  
RUC83



13 mm from pivot to tip, left, 105 mm long.




**General Scissors Sharp Pointed**

**Code**  
RUG01



Sharp pointed 9.5 cm scissors.




**General Scissors Blunt End**

**Code**  
RUG06



Blunt Ended 9 cm scissors.




**General Scissors Long Blunt End**

**Code**  
RUG14




Blunt ended 13 cm scissors.



**Steven Scissors Curved**

**Code**  
RUG23





A series of 20 horizontal blue lines spanning the width of the page, providing a space for taking notes.



---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---






---

---


---

---


---

<b>Forceps</b>			
	Castroviejo Tying Forceps	1401031 / 1401032	1
	Curved Tying Forceps	1401018 / 1401026	1
	Girard Tying Forceps	1401041 / 1401042	1
	McPherson Tying Forceps	1401013 / 1401014 / 1401015 / 1401016	1
	McPherson Tying Forceps	1401023 / 1401024	1
	Straight Tying Forceps	1401011 / 1401012	2
	Straight Tying Forceps	1401021 / 1401022	2
	Bishop-Harmon Tissue Forceps	1401156	2
	Bonn Toothed Forceps	1401160	2
	Colibri Conjunctiva Toothed Forceps	1401117	2
	Colibri Toothed Forceps	1401113 / 1401123	3
	Colibri Toothed Forceps	1401153	3
	Harms-Colibri Toothed Forceps	1401151 / 1401161	3
	Kelman-McPherson Toothed Forceps	1401114	3
	Polack Corneal Suturing Forceps	1401372	3
	Polack Corneal Suturing Forceps	1401381	4
	Straight Toothed Forceps	1401111 / 1401112	4
	Straight Toothed Forceps	1401121	4
	Straight Toothed Forceps	1401122 / 1401141	4
	Curved Toothed Forceps	1401152 / 1401154	4
	Tissue Forceps	1401159 / 1401162	5
	Toothed Forceps	1401131 / 1401132	5
	Angled Notched Forceps	1401063	5
	Curved Notched Forceps	1401065	5
	Hofmann-Polack Corneal Suturing Forceps	1401371	5
	Pierse Corneal Forceps	1401081 / 1401082	6
	Pierse-Colibri Corneal Forceps	1401101	6
	Straight Notched Forceps	1401061 / 1401062	6
	Straight Notched Forceps	1401071 / 1401072	6
	Bloomberg Pierse Capsulorhexis Forceps	1401222	7
	Capsulorhexis Forceps	1401225	7
	Capsulorhexis Forceps	1401226	7
	Corydon Capsulorhexis Forceps	1401223	7
	Corydon Style Capsulorhexis Forceps	1401176 / 1401186	7
	Lehner-Utrata Style Capsulorhexis Forceps	1401221	8
	O'Gawa-Utrata Capsulorhexis Forceps	1401224	8
	Utrata Style Capsulorhexis Forceps	1401173 / 1401174	8
	Utrata Style Capsulorhexis Forceps	1401227	8
	Lens Forceps	1401261 / 1401271	9
	Kelman-McPherson Lens Forceps	1401301	9
	Kratz Style Lens Forceps	1401283	9
	Chalazion Forceps	1401405 / 1401406 / 1401407 / 1401408	9
	Baby Jones Towel Clamps	1401417	10
	Cilia Forceps	1401412	10
	Hartmann Curved Hemostatic Mosquito Forceps	1401415	10
	Hartmann Straight Hemostatic Mosquito Forceps	1401414	10
	Scleral Plug Forceps	1401416	10
	Jeweller Style Forceps	1401421	11
	Angled Jeweller Style Forceps	1401422	11
	Angled Jeweller Style Forceps (45 degree)	1401423	11
	Curved Jeweller Style Forceps	1401424	11
	<b>Speculae</b>		
	Barraquer Solid Speculum	1405008	12
	Eye Speculum	1405014	12
	Kratz Barraquer Titanium Wire Speculum	1405007	12
	Lieberman Speculum, K-wire	1405005	12
	Lieberman Speculum, V-wire	1405002	12
	Naqahara Titanium Wire Speculum	1405009	13
	Reversible Lid Speculum	1405006	13
	Shapira Speculum	1405012	13
	Slade LASIK Titanium Wire Speculum	1405003	13
	Titanium Wire Speculum	1405001	13
	Titanium Wire Speculum	1405015	14
	<b>Nucleus Rotators &amp; Choppers</b>		
	Akahoshi Nucleus Stainer/manipulator	1404140	15
	Claw Phaco Chopper	1404130	15
	Drysdale Nucleus Manipulator	1404113	15
	Haefliger Phaco Cleaver	1404139	15
	Helveston Retractors	1404158 / 1404159 / 1404160	15
	Helveston Scleral Marking Ruler	1404161	16
	Intracapsular Manipulator	1404120	16
	Intracapsular Manipulator	1404151	16
	Jaffe-Bechert Lens Nucleus Rotator	1404105	16
	Lester IOL Manipulator	1404106 / 1404107	16
	Manipulator	1404135	17
	Nichamin Phaco Chopper	1404123 / 1404124	17
	Phaco Chopper	1404121 / 1404122	17
	Phaco Chopper	1404125 / 1404125L	17
	Phaco Chopper	1404126 / 1404126L	17
	Phaco Chopper	1404136 / 1404137	18
	Phaco Chopper	1404138 / 1404138L	18
	Phaco Chopper	1404141	18
	Phaco Chopper	1404142 / 1404142L	18
	Phaco Chopper	1404302	18
	Phaco Chopper	1404303	19
	Repositor	1404114 / 1404115	19
	Schocket Double Ended Scleral Depressor	1404304	19
	<b>Hooks</b>		
	Helveston Finder Hook	1404155 / 1404156	20
	Helveston Muscle Hook	1404152 / 1404153 / 1404154	20
	Helveston Teaser Hook	1404157	20
	Kuglen Iris Hook	1404101 / 1404102	20
	Retinal Detachment Hook	1404119	20
	Sinsky Lens Manipulating Hook	1404110	21
	<b>Needle Holders</b>		
	Barraquer Needle Holder	1402238 / 1402245	22
	Barraquer Needle Holder-Curved	1402101 / 1402106	22
	Bechart-Sinsky Needle Holder	1402235 / 1402243	22
	Castroviejo Needle Holder	1402191 / 1402192	22
	Castroviejo Needle Holder	1402193 / 1402195	22
	Castroviejo Needle Holder	1402204 / 1402205	23
	Castroviejo Needle Holder	1402203 / 1402206	23
	Castroviejo Needle Holder	1402211 / 1402212	23
	Castroviejo Needle Holder	1402213 / 1402214	23
	Curved Needle Holder	1402011 / 1402016	23
	Curved Needle Holder	1402233 / 1402241	24
	Fine Curved Needle Holder	1402237 / 1402244	24
	Side Curved Needle Holder	1402236	24
	Straight Needle Holder	1402234 / 1402242	24
	Troutman Needle Holder-Curved	1402231 / 1402232	24


### Fixation Instruments

	Fine Fixation Ring <a href="#">1404149</a>	25
	Fine/Thornton Fixation Ring <a href="#">1404147 / 1404148</a>	25
	Fixation Forceps <a href="#">1401373 / 1401375</a>	25
	Helveston Fixation Forceps <a href="#">1401382</a>	25
	Hofman-Thornton Fixation Ring <a href="#">1404146</a>	25


### Lasik

	Buratto Flap Protector <a href="#">1405031</a>	26
	Corneal Epithelium Spatula <a href="#">1405027</a>	26
	Flap Elevator <a href="#">1405022</a>	26
	Genisi Speculum LASIK <a href="#">1405034</a>	26
	Katzen LASIK Flap Unzipper <a href="#">1405046</a>	26
	LASIK Flap Manipulator <a href="#">1405045</a>	27
	Lieberman Solid Speculum <a href="#">1405029</a>	27
	Machat Adjustable Speculum <a href="#">1405030</a>	27
	Steinert LASIK Speculum LASIK <a href="#">1405033</a>	27

### Other

	Blade Breaker <a href="#">1410001</a>	28
	Blade Handle <a href="#">1410002</a>	28
	Braunstein Fixed Caliper <a href="#">1410003 / 1410004</a>	28
	Castroviejo Caliper <a href="#">1410005</a>	28
	Christen Punctum Punch <a href="#">1410014</a>	28
	Descemet Membrane Punch <a href="#">1410012</a>	29
	Holth Corneo scleral Punch <a href="#">1410013</a>	29
	Irrigating Vectis (3x8) <a href="#">1410302</a>	29
	Luntz-Dodick Trabulectomy Punch <a href="#">1410011 / 1410015 / 1410016</a>	29
	Trabeculotome set <a href="#">1404134</a>	29
	Trabeculotome set <a href="#">1404134C</a>	30
	Urrets-Zavalía sclera marker <a href="#">1404162</a>	30

### Stainless Steel Scissors

	Westcott Scissors <a href="#">UA11</a>	31
	Ong Type Scissors <a href="#">RUC80</a>	31
	Vannas Scissors <a href="#">RUC81</a>	31
	Castroviejo Corneal Section Scissors (Right) <a href="#">RUC82</a>	31
	Castroviejo Corneal Section Scissors (Left) <a href="#">RUC83</a>	31
	General Scissors Sharp Pointed <a href="#">RUG01</a>	32
	General Scissors Blunt End <a href="#">RUG06</a>	32
	General Scissors Long Blunt End <a href="#">RUG14</a>	32
	Steven Scissors Curved <a href="#">RUG23</a>	32



**Distributed by:**

**Medix21 Surgical**

29 Jarden Mile Ngauranga 6035

p 0800 633 4921 / 021 734 762

e [camille@medix21.co.nz](mailto:camille@medix21.co.nz)

w [medix21surgical.co.nz](http://medix21surgical.co.nz)

**Surgitrac Instruments Limited**

1213 Wing On Plaza, 62 Mody Road, Tsim Sha Tsui East, Kowloon, Hong Kong.

Tel: +44(0)161 776 7626 Fax: +44(0)161 776 7616 Email: [info@surgitrac.com](mailto:info@surgitrac.com) [www.surgitrac.com](http://www.surgitrac.com)