





***auro**Trocar*

Self Sealing Ophthalmic Trocar (Valved Device)



Aurotrocar, is a self sealing ophthalmic trocar, designed and manufactured by Aurolab, a division of Aravind Eye Care System—a name synonymous with quality and commitment. Developed as an affordable product, to reduce the burden of high costs to both hospitals and patients, this valved device boasts of excellent features that meet the needs of retinal surgeons.



FEATURES

- Sharp needle for smooth and high precision incision
- Tapered point cannula for smooth penetration
- Cone shaped cannula for easy insertion of instruments
- Self sealing silicone valve cap prevents leaking
- Low friction silicone valves designed for smooth instrument exchange
- Precise scleral marking point for safe incision positioning in 3 mm & 4 mm
- Perfect locking of adaptor with the cannula to avoid choroidal damage
- Color coded handles for 23GA and 25GA
- Double packed in blister and pouch

OTHER INFORMATION

- One box contains – 3 valved trocar cannulas mounted on sharp blades with one infusion line
- One step trocar
- Sizes – 23GA & 25GA



1, Sivagangai Main Road, Veerapanjan, Madurai 625 020, India
Phone: 91 452 3096100 · Fax: 91 452 2446200 · E-mail: info@aurolab.com · Web: www.aurolab.com